


Product Datasheet

CSDE1 antibody (orb373315)

Description	Rabbit polyclonal antibody to CSDE1
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IP, WB
Immunogen	Recombinant fusion protein of human CSDE1(NP_009089.4).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	89KDa
Uniprot ID	075534
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...