



Product Datasheet

ELAVL3 antibody (orb373336)


Description	Rabbit polyclonal antibody to ELAVL3
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide of human ELAVL3(NP_001411.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	37KDaa
Uniprot ID	Q14576
Dilution Range	WB: 1:500-1000, IHC-P: 1:50-200, IP: 1:20-50R, IP: 1:20-50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemical staining of rat brai...



Immunohistochemical staining of mouse br...