



Product Datasheet

KHDRBS2 antibody (orb373338)


Description	Rabbit polyclonal antibody to KHDRBS2
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide of human KHDRBS2(NP_689901.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	39kDa
Uniprot ID	Q5VWX1
Dilution Range	WB: 1:500-1000, IHC-P: 1:100-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemical staining of human pr...



Immunohistochemical staining of human es...