


## Product Datasheet

### TNRC6A antibody (orb373341)

|                            |                                                    |
|----------------------------|----------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to TNRC6A               |
| <b>Reactivity</b>          | Human, Mouse, Rat                                  |
| <b>Conjugation</b>         | Unconjugated                                       |
| <b>Tested Applications</b> | ICC, IF, WB                                        |
| <b>Immunogen</b>           | A synthetic peptide of human TNRC6A(NP_055309.2).  |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH 7.3. |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.        |
| <b>Note</b>                | For research use only                              |
| <b>Isotype</b>             | IgG                                                |
| <b>Clonality</b>           | Polyclonal                                         |
| <b>MW</b>                  | 180kDa                                             |
| <b>Uniprot ID</b>          | <b>Q8NDV7</b>                                      |
| <b>Dilution Range</b>      | WB: 1:500-1000                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                    |



Western blot analysis of extracts of var...