



Product Datasheet

LMX1B antibody (orb373377)

Description	Rabbit polyclonal antibody to LMX1B
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human LMX1B(NP_002307.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	53KDa
Uniprot ID	O60663
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of BT-474 cells.



Immunohistochemical staining of human liver tissue.