



## Product Datasheet

### PIP antibody (orb373380)


|                            |                                                       |
|----------------------------|-------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to PIP                     |
| <b>Reactivity</b>          | Human, Mouse, Rat                                     |
| <b>Conjugation</b>         | Unconjugated                                          |
| <b>Tested Applications</b> | ICC, IF, IHC                                          |
| <b>Immunogen</b>           | Recombinant fusion protein of human PIP(NP_002643.1). |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.    |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.           |
| <b>Note</b>                | For research use only                                 |
| <b>Isotype</b>             | IgG                                                   |
| <b>Clonality</b>           | Polyclonal                                            |
| <b>MW</b>                  | 16kDa                                                 |
| <b>Uniprot ID</b>          | <b>P12273</b>                                         |
| <b>Dilution Range</b>      | WB: 1:500-2000, IHC-P: 1:50-100, IF/ICC: 1:50-100     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                       |



Immunohistochemical staining of rat lung...



Immunohistochemical staining of human pr...



Immunohistochemical staining of mouse lu...