



Product Datasheet

TSPAN7 antibody (orb373388)


Description	Rabbit polyclonal antibody to TSPAN7
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human TSPAN7(NP_004606.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	28KDa
Uniprot ID	P41732
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of Mou...



Immunofluorescence analysis of Mouse pan...



Immunofluorescence analysis of Rat Pancr...