



Product Datasheet

AGPAT2 antibody (orb373410)

Description	Rabbit polyclonal antibody to AGPAT2
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human AGPAT2(NP_006403.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	31kDa
Uniprot ID	O15120
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of HT-29 cells.



Immunofluorescence analysis of HeLa cells.