



Product Datasheet

AKAP5 antibody (orb373411)


Description	Rabbit polyclonal antibody to AKAP5
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human AKAP5(NP_004848.3).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	47kDa
Uniprot ID	P24588
Dilution Range	WB: 1:500-2000, IHC-P: 1:100-200, IF/ICC: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse brain tissue.



Immunohistochemical staining of rat brain tissue.



Immunofluorescence analysis of HeLa cells.