


Product Datasheet

CLDN2 antibody (orb373422)

Description	Rabbit polyclonal antibody to CLDN2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	Recombinant fusion protein of human CLDN2 (NP_001164566.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P57739
Entrez	9075
Dilution Range	IF/ICC: 1:50-100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...