



Product Datasheet

HPS1 antibody (orb373434)


Description	Rabbit polyclonal antibody to HPS1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Recombinant fusion protein of human HPS1(NP_872577.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	36kDa/75kDa/76kDa/79kDa
Uniprot ID	Q92902
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100, IF/ICC: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse st...



Immunohistochemical staining of rat kidn...



Immunohistochemical staining of human st...