


Product Datasheet

SETMAR antibody (orb373465)

Description	Rabbit polyclonal antibody to SETMAR
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human SETMAR(NP_006506.3).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	40kDa/62kDa/78kDa
Uniprot ID	Q53H47
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:20-50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...