



Product Datasheet

WFDC2 antibody (orb373476)


Description	Rabbit polyclonal antibody to WFDC2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human HE4/WFDC2(NP_006094.3).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	20KDa
Uniprot ID	Q14508
Dilution Range	IHC-P: 1:100-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat hear...



Immunohistochemical staining of human pr...



Immunohistochemical staining of human st...