



Product Datasheet

TARS antibody (orb373519)


Description	Rabbit polyclonal antibody to TARS
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human TARS(NP_689508.3).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	83kDa
Uniprot ID	P26639
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100, IF/ICC: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse ki...



Immunohistochemical staining of rat hear...



Immunohistochemical staining of mouse br...