



Product Datasheet

VSNL1 antibody (orb373520)

Description	Rabbit polyclonal antibody to VSNL1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human VSNL1/VILIP-1(NP_003376.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	18KDa
Uniprot ID	P62760
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse brain...



Western blot analysis of extracts of var...