


Product Datasheet

SH3GLB1 antibody (orb373536)

Description	Rabbit polyclonal antibody to SH3GLB1
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human Bif-1(NP_057093.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	29kDa/40kDa/43kDa
Uniprot ID	Q9Y371
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of A549 cell...