



Product Datasheet

TRMT1 antibody (orb373540)

Description	Rabbit polyclonal antibody to TRMT1
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human TRMT1(NP_060192.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	90kDa
Uniprot ID	Q9NXH9
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human lung tissue.



Immunohistochemical staining of human gastric tissue.