



Product Datasheet

ALKBH8 antibody (orb373549)


Description	Rabbit polyclonal antibody to ALKBH8
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human ALKBH8 (NP_620130.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	75kDa
Uniprot ID	Q96BT7
Entrez	91801
Dilution Range	WB: 1:500-2000, IHC-P: 1:100-200, IF/ICC: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse lung tissue.



Immunohistochemical staining of human gastric tissue.



Immunofluorescence analysis of A549 cells.