



Product Datasheet

NCR3 antibody (orb373553)


Description	Rabbit polyclonal antibody to NCR3
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human NCR3(NP_667341.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	30kDa
Uniprot ID	O14931
Dilution Range	WB: 1:500-2000, IHC-P: 1:100-200, IF/ICC: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse lung tissue.



Immunohistochemical staining of rat lung tissue.



Immunohistochemical staining of human esophagus tissue.