


Product Datasheet

SLC35A2 antibody (orb373579)

Description	Rabbit polyclonal antibody to SLC35A2
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Recombinant fusion protein of human SLC35A2(NP_001027460.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	25kDa/34kDa/41kDa/44kDa
Uniprot ID	P78381
Dilution Range	IF :50-200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...