


Product Datasheet

Acetyl-Histone H3-K27 antibody (orb373581)

Description	Rabbit polyclonal antibody to Acetyl-Histone H3-K27
Reactivity	All, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, ICC, IF, IHC, IP, WB
Immunogen	A synthetic acetylated peptide around K27 of human Histone H3(NP_003520.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	17KDa
Uniprot ID	P68431
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200, IP: 1:50-200, ChIP: 1:20-100, ChIP
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...