

Product Datasheet

Acetyl-Histone H4-K8 antibody (orb373586)

Description Rabbit polyclonal antibody to Acetyl-Histone H4-K8

Reactivity All, Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ChIP, ICC, IF, IHC, IP, WB

Immunogen A synthetic acetylated peptide around K8 of human Histone H4(NP_003529.1).

Preservatives PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Storage Store at -20°C. Avoid freeze / thaw cycles.

Note For research use only


Isotype IgG

MW 11kDa


Uniprot ID **P62805**

Dilution Range WB: 1:1000-3000, IHC-P: 1:200-500, IF/ICC: 1:500-1000, IP: 1:200-500, ChIP: 1:20-50


Expiration Date 12 months from date of receipt.



Western blot analysis of extracts of var...



Immunoprecipitation analysis extracts of...



Immunofluorescence analysis of HeLa cell...