



Product Datasheet

PNLIP antibody (orb373629)


Description	Rabbit polyclonal antibody to PNLIP
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human PNLIP(NP_000927.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	51kDa
Uniprot ID	P16233
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of BxPC3 cell lysa...



Immunohistochemical staining of human li...



Immunofluorescence analysis of U2OS cell...