


Product Datasheet

CYB5R3 antibody (orb373633)

Description	Rabbit polyclonal antibody to CYB5R3
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human CYB5R3(NP_001165131.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	34-36kDa
Uniprot ID	P00387
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
extracts of
var...