



Product Datasheet

TTBK2 antibody (orb373650)

Description	Rabbit polyclonal antibody to TTBK2
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human TTBK2 (NP_775771.3).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	137kDa
Uniprot ID	Q6IQ55
Entrez	146057
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:20-50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of rat...



Immunohistochemical staining of human ga...