


Product Datasheet

MLXIPL antibody (orb373656)

Description	Rabbit polyclonal antibody to MLXIPL
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human MLXIPL(NP_116569.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	72kDa\110kDa
Uniprot ID	Q9NP71
Dilution Range	WB: 1:500-2000, IP: 1:10-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...