



## Product Datasheet

### SLC2A4 antibody (orb373661)

|                            |                                                  |
|----------------------------|--------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to SLC2A4             |
| <b>Reactivity</b>          | Human, Mouse, Rat                                |
| <b>Conjugation</b>         | Unconjugated                                     |
| <b>Tested Applications</b> | ICC, IF, IHC, WB                                 |
| <b>Immunogen</b>           | A synthetic peptide of human GLUT4(NP_001033.1). |
| <b>Preservatives</b>       | PBS with 0.05% proclin300,50% glycerol,pH7.3.    |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.      |
| <b>Note</b>                | For research use only                            |
| <b>Isotype</b>             | IgG                                              |
| <b>Clonality</b>           | Polyclonal                                       |
| <b>MW</b>                  | 55KDs                                            |
| <b>Uniprot ID</b>          | <b>P14672</b>                                    |
| <b>Dilution Range</b>      | WB: 1:500-2000, IHC-P: 1:50-200                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                  |



Western blot analysis of extracts of var...



Immunohistochemical staining of mouse he...