


Product Datasheet

CTSV antibody (orb373678)

Description	Rabbit polyclonal antibody to CTSV
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human CTSV (NP_001324.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	37kDa
Uniprot ID	O60911
Entrez	1515
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...