



Product Datasheet

MAP2K4 antibody (orb373717)

Description	Rabbit polyclonal antibody to MAP2K4
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic peptide of human MAP2K4(NP_003001.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	44KDa
Uniprot ID	P45985
Dilution Range	WB: 1:1000-4000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemical staining of human ki...