



Product Datasheet

PDE4DIP antibody (orb373749)


Description	Rabbit polyclonal antibody to PDE4DIP
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human PDE4DIP(NP_071754.3).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	265kDa
Uniprot ID	Q5VU43
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of mouse skeletal muscle.



Immunohistochemical staining of rat heart tissue.



Immunohistochemical staining of human kidney tissue.