



Product Datasheet

COQ7 antibody (orb373751)


Description	Rabbit polyclonal antibody to COQ7
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human COQ7(NP_057222.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	26kDa
Uniprot ID	Q99807
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of mouse skeletal muscle and brain.



Immunohistochemical staining of mouse liver tissue.



Immunohistochemical staining of human liver tissue.