



Product Datasheet

AARS2 antibody (orb373791)

Description	Rabbit polyclonal antibody to AARS2
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human AARS2(NP_065796.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	107kDa
Uniprot ID	Q5JZ9
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of various cell lines.



Immunohistochemical staining of human liver tissue.