


Product Datasheet

KRT10 antibody (orb373847)

Description	Rabbit polyclonal antibody to KRT10
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human KRT10(NP_000412.4).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	60KDa
Uniprot ID	P13645
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...