



Product Datasheet

ASCC3 antibody (orb373883)

Description	Rabbit polyclonal antibody to ASCC3
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Recombinant fusion protein of human ASCC3(NP_071374.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	13kDa/83kDa/251kDa
Uniprot ID	Q8N3C0
Dilution Range	WB: 1:500-2000, IF/ICC: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunofluorescence analysis of A549 cell...