


## Product Datasheet

### LAMB3 antibody (orb373966)

|                            |                                                           |
|----------------------------|-----------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to LAMB3                       |
| <b>Species/Host</b>        | Rabbit                                                    |
| <b>Reactivity</b>          | Human                                                     |
| <b>Conjugation</b>         | Unconjugated                                              |
| <b>Tested Applications</b> | IHC                                                       |
| <b>Immunogen</b>           | Recombinant fusion protein of human LAMB3 (NP_000219.2).  |
| <b>Preservatives</b>       | Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3. |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.               |
| <b>Note</b>                | For research use only                                     |
| <b>Isotype</b>             | IgG                                                       |
| <b>Clonality</b>           | Polyclonal                                                |
| <b>Purity</b>              | Affinity purification                                     |
| <b>Uniprot ID</b>          | <b>Q13751</b>                                             |
| <b>Entrez</b>              | <b>3914</b>                                               |
| <b>Dilution Range</b>      | WB: 1:500-2000, IHC-P: 1:50-100                           |
| <b>Expiration Date</b>     | 12 months from date of receipt.                           |



Immunohistochemical staining of human st...