



Product Datasheet

PLCB2 antibody (orb373973)

Description	Rabbit polyclonal antibody to PLCB2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human PLCB2(NP_004564.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	140KDa
Uniprot ID	Q00722
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extract of vari...



Immunohistochemical staining of human st...