


Product Datasheet

STXBP3 antibody (orb373979)

Description	Rabbit polyclonal antibody to STXBP3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human STXBP3 (NP_009200.2).
Preservatives	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	68kDa
Uniprot ID	O00186
Entrez	6814
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
extract of
vari...