


Product Datasheet

ADAM33 antibody (orb374022)

Description	Rabbit polyclonal antibody to ADAM33
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human ADAM33(NP_001269376.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	87kDa
Uniprot ID	Q9BZ11
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
extracts of
var...