



Product Datasheet

ATAD3B antibody (orb374025)


Description	Rabbit polyclonal antibody to ATAD3B
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human ATAD3B(NP_114127.3).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	68KDa
Uniprot ID	Q5T9A4
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extract of vari...



Immunohistochemical staining of mouse heart...



Immunohistochemical staining of human li...