



Product Datasheet

PPP4R1 antibody (orb374074)


Description	Rabbit polyclonal antibody to PPP4R1
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Recombinant fusion protein of human PPP4R1(NP_001035847.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	130kDa
Uniprot ID	Q8TF05
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemical staining of human li...



Immunohistochemical staining of human st...