



Product Datasheet

TMEM189 antibody (orb374090)


Description	Rabbit polyclonal antibody to TMEM189
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide of human TMEM189 (NP_954580.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	31kDa
Uniprot ID	A5PLL7
Entrez	387521
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of MCF...



Immunohistochemical staining of human li...



Immunohistochemical staining of human st...