



Product Datasheet

YTHDF3 antibody (orb374095)


Description	Rabbit polyclonal antibody to YTHDF3
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IP, WB
Immunogen	Recombinant fusion protein of human YTHDF3(NP_689971.4).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	64KDa
Uniprot ID	Q7Z739
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human esophagus tissue.



Immunohistochemical staining of rat brain tissue.



Immunohistochemical staining of mouse stomach tissue.