



Product Datasheet

HRG antibody (orb374122)


Description	Rabbit polyclonal antibody to HRG
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human HRG(NP_000403.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	80kDa
Uniprot ID	P04196
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of mouse pancreas.



Immunohistochemical staining of rat lung tissue.



Immunohistochemical staining of mouse lung tissue.