



Product Datasheet

MTERFD3 antibody (orb374181)


Description	Rabbit polyclonal antibody to MTERFD3
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human MTERFD3(NP_079474.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	44KDa
Uniprot ID	Q49AM1
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemical staining of mouse ki...



Immunohistochemical staining of rat kidn...