



## Product Datasheet

### KCNK9 antibody (orb374220)

|                            |                                                     |
|----------------------------|-----------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to KCNK9                 |
| <b>Reactivity</b>          | Human, Mouse, Rat                                   |
| <b>Conjugation</b>         | Unconjugated                                        |
| <b>Tested Applications</b> | WB                                                  |
| <b>Immunogen</b>           | A synthetic peptide of human KCNK9(NP_001269463.1). |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.  |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.         |
| <b>Note</b>                | For research use only                               |
| <b>Isotype</b>             | IgG                                                 |
| <b>Clonality</b>           | Polyclonal                                          |
| <b>MW</b>                  | 42kDa/53kDa                                         |
| <b>Uniprot ID</b>          | <b>Q9NPC2</b>                                       |
| <b>Dilution Range</b>      | WB: 1:500-2000                                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                     |



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...