


Product Datasheet

ALDH7A1 antibody (orb374221)

Description	Rabbit polyclonal antibody to ALDH7A1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human ALDH7A1(NP_001173.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	58kDa
Uniprot ID	P49419
Dilution Range	WB: 1:1000-4000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
extracts of
var...