


Product Datasheet

GTF2E2 antibody (orb374223)

Description	Rabbit polyclonal antibody to GTF2E2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human GTF2E2(NP_002086.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	33KDa
Uniprot ID	P29084
Dilution Range	WB: 1:200-2000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
extracts of
var...