


Product Datasheet

NUMB antibody (orb374235)

Description	Rabbit polyclonal antibody to NUMB
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human NUMB(NP_003735.3).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	71KDa
Uniprot ID	P49757
Dilution Range	WB: 1:500-2000, IHC-P: 1:20-200, IF/ICC: 1:10-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...