



Product Datasheet

LCN1 antibody (orb374257)

Description	Rabbit polyclonal antibody to LCN1
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human LCN1(NP_002288.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	16kDa
Uniprot ID	P31025
Dilution Range	WB: 1:200-500, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemical staining of human pr...