


## Product Datasheet

### ANKH antibody (orb374286)

|                            |                                                        |
|----------------------------|--------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to ANKH                     |
| <b>Reactivity</b>          | Human                                                  |
| <b>Conjugation</b>         | Unconjugated                                           |
| <b>Tested Applications</b> | WB                                                     |
| <b>Immunogen</b>           | Recombinant fusion protein of human ANKH(NP_473368.1). |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide,50% glycerol,pH7.3.        |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.            |
| <b>Note</b>                | For research use only                                  |
| <b>Isotype</b>             | IgG                                                    |
| <b>Clonality</b>           | Polyclonal                                             |
| <b>MW</b>                  | 54kDa                                                  |
| <b>Uniprot ID</b>          | <b>Q9HCJ1</b>                                          |
| <b>Dilution Range</b>      | WB: 1:500-2000, IHC-P: 1:50-100                        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                        |



Immunohistochemical staining of human pr...